

12th Webinar

Do you know the high sensitive ToxProtect?



Fish in a tank - best choice for a simple, affordable and effective Early Warning System

Welcome



Ylva Tischler

Detlev Lohse

bbe Team

12th Webinar

Part 1 **Why do we measure Fish toxicity?**

Part 2 **What is the Value of *ToxProtect*?**

Part 3 **Construction & Function**



Part 4 **Application**

Conclusion

Feedback

Follow up

The Testorganism in the Phylogenetic Tree



OPPTS 850.1075 EPA

Test: ACUTE TOXICITY FOR FISH

The purpose of this test is to determine the acute lethal toxicity of a substance to fish in fresh water and marine water.



Literature on Fish Toxicity

Welcome to Search and Browse for EXTOXNET...

Pesticide Information Profiles (PIPs)

EXTOXNET is a cooperative effort of University of California-Davis, Oregon State University, Michigan State University, Cornell University, and the University of Idaho.
Primary files are maintained and archived at Oregon State University.


Browse... Browse these PIPs by 'clicking' on the title:

- [WHAT ARE PIPs?](#)
- [2,4-D - EXTOXNET PIP](#)
- [2,4-DB - EXTOXNET PIP](#)
- [4-AMINOPYRIDINE - EXTOXNET PIP](#)
- [DCPA, CHLORTHAL, CLORTHAL-DIMETHYL - EXTOXNET PIP](#)
- [MCPA - EXTOXNET PIP](#)
- [ABAMECTIN - EXTOXNET PIP](#)
- [ACEPHATE \(ORTHENE\) - EXTOXNET PIP](#)
- [ACETOCHLOR - EXTOXNET PIP](#)
- [ACIFLUORFEN - EXTOXNET PIP](#)
- [ALACHLOR - EXTOXNET PIP](#)
- [ALDICARB - EXTOXNET PIP](#)
- [ALLETHRIN - EXTOXNET PIP](#)
- [ALUMINUM PHOSPHIDE - EXTOXNET PIP](#)
- [AMETRYN - EXTOXNET PIP](#)
- [AMITRAZ - EXTOXNET PIP](#)
- [AMITROLE - EXTOXNET PIP](#)
- [AMMONIUM SULFAMATE - EXTOXNET PIP](#)
- [ATRAZINE - EXTOXNET PIP](#)


.....may cause fish and aquatic invertebrate deaths [43].

Chlorpyrifos toxicity to fish may be related to water temperature. The 96-hour LC50 for chlorpyrifos is 0.009 mg/L in mature rainbow trout, 0.098 mg/L in lake trout, 0.806 mg/L in goldfish, 0.01 mg/L in bluegill, and 0.331 mg/L in fathead minnow [50].....

How does the Static Test work?



Tube #	1	2	3	4	5
volume of water	0	1 mL	1 mL	1 mL	1 mL
dilution ratio	1	1/2	1/4	1/8	1/16



Is the Static Test an Option?

No




Why not?



Online Toxicity and Pollution Wave

Toximeter signal —
 Pollutant concentration - - -



The Alternative – Chemical Analysis?



ToxProtect 64

ToxProtect is a cost-effective, automated biomonitor for rapid detection of acute toxic substances in water

- detection of fish movements by an array of light barriers
- remote access to data including all optional sensors and online operation
- low maintenance
- simple installation




What makes a *ToxProtect* Attractive in Toxicity Monitoring?

- **Continuous monitoring**
- **Ease of operation**
- **Sensitive vertebrate – related physiology**
- **Reliability of hardware**
- **Wide temperature range: cold & warm water fish**



Set-up of the *ToxProtect*



Essential Parts of the *ToxProtect*



Essential Parts of the *ToxProtect*



To put the Items together

■ **Hardware ToxProtect**



■ **Fish/Testorganism**



■ **AlarmSoftware**

To each Place the Suitable Fish



Medaka



Zebrafish



Dace (Leuciscus)



Tigerbarb



Blue gill



Fathead minnow

Selection for the Type of Substances

www.pesticideinfo.org/List_EcoChemSpecies.jsp?Taxa_Group='Fish'

PAN Pesticides Database - Chemical Toxicity Studies on Aquatic Organisms

[Home](#) > [Ecotoxicity Search](#)

Aquire Search Results

The 'Fish' organism group contains the following species. Click on any species listed below for a list of studies. Note that if we do not find any results in the results. The PAN Pesticide database only includes pesticide-related chemicals.



Common Name	Scientific Name	Number of Studies
Blue Bream	Abramis ballerus	4
Bream	Abramis brama	28
Yellowfin goby	Acanthogobius flavimanus	5
Surf bream	Acanthopagrus australis	3
Porgy	Acanthopagrus schlegelii	36
Minnow, tanago	Achellognathus morioakae	15
Siberian Sturgeon	Acipenser baerii	47
Lake sturgeon	Acipenser fulvescens	38
Sturgeon	Acipenser gueldenstaedti colch	1
Russian sturgeon	Acipenser gueldenstaedti	2
Sterlet	Acipenser ruthenus	1
Sevruga, stellate sturgeon	Acipenser stellatus	8
White sturgeon	Acipenser transmontanus	42
Sturgeon family	Acipenseridae	1
Minnow	Acrossochellus paradoxus	13
Spiny rayed fish class	Actinopterygii	3
Diamond killifish	Adinia xenica	3
Hooknose	Agonus cataphractus	48
Longfin dace	Agosia chrysoaster	3
Bleak	Alburnus albidus	5
Bleak	Alburnus alburnus	119
Yelloweye mullet	Aldrichetta forsteri	9
Blueback herring	Alosa aestivalis	9
Alewife	Alosa pseudoharengus	5
American shad	Alosa sapidissima	7
Glassy, Perchlet	Ambassis commersoni	3
Bald glassy	Ambassis gymnocephalus	3

(Carboxymethoxy)butanedioic acid, Trisodium salt [Show \(Carboxymethoxy\)butanedioic acid, Trisodium salt studies for all species](#)

Species	Endpoint	Effect	Substance	Duration	LC50	LC10	LC5	NOEC	Concentration	Route	Other	Year	Reference	
Zebra danio Danio rerio	Mortality	Mortality	FRY, 3 WK	96 h	LC50	2,100,000	1,700,000	2,500,000	ug/L	A	63.9 % AI	Static	Not Acutely Toxic	1982 J. Environ. Qual. of Life Report No. EUR 7549:284-295
Zebra danio Danio rerio	Mortality	Mortality	JUVENILE, 12 WK	96 h	LC50	2,600,000	2,500,000	2,700,000	ug/L	A	63.9 % AI	Static	Not Acutely Toxic	1982 J. Environ. Qual. of Life Report No. EUR 7549:284-295
Zebra danio Danio rerio	Mortality	Mortality	ADULT, 20 WK	96 h	LC50	2,500,000	1,800,000	3,100,000	ug/L	A	63.9 % AI	Static	Not Acutely Toxic	1982 J. Environ. Qual. of Life Report No. EUR 7549:284-295


(T-4)-bis(2-Hydroxypropanoato-0'1, O2)zinc [Show \(T-4\)-bis\(2-Hydroxypropanoato-0'1, O2\)zinc studies for all species](#)

Species	Endpoint	Effect	Substance	Duration	LC50	LC10	LC5	NOEC	Concentration	Route	Other	Year	Reference	
Zebra danio Danio rerio	Mortality	Mortality	NR	96 h	LC50	23,100	-	-	ug/L	T	NR	Static	Slightly Toxic	1998 Chemosphere 37(7):1317-1333
Zebra danio Danio rerio	Mortality	Mortality	NR	96 h	NOEC	12,900	-	-	ug/L	T	NR	Static		1998 Chemosphere 37(7):1317-1333

1,1,1-trichloroethane [Show 1,1,1-trichloroethane studies for all species](#)


Species	Endpoint	Effect	Substance	Duration	LC50	LC10	LC5	NOEC	Concentration	Route	Other	Year	Reference	
Zebra danio Danio rerio	Behavior	Behavioral changes, general	NR	14 d	NOEC	3,400	-	-	ug/L	F	NR	Flow through	1990 Testbericht: Wassergefährdende Stoffe Fraunhofer-Institut für Umweltchemie und Ökotoxikologie, Schmallenberg (OECD Data File)	
Zebra danio Danio rerio	Mortality	Mortality	NR	48 h	LC50	79,000	-	-	ug/L	F	NR	Flow through	Slightly Toxic	1990 Testbericht: Wassergefährdende Stoffe Fraunhofer-Institut für Umweltchemie und Ökotoxikologie, Schmallenberg (OECD Data File)

Response Time



To put the Items together

- **Hardware ToxProtect**
- **Fish/Testorganism**
- **Alarm Software**




Alarm Recognition

Activity

Specific activity

Thresholds

Alarm



False Alarm Prevention

Problem:

- false alarms occur due to the natural variation in behaviour
- due to the low risk of an event, alarms can be wrong

Consequence:



- an alarm evaluation is required to reduce the risk

Solutions

- second independent instrument
- second type of test in the same instrument



Alarm Verification = False Alarm Prevention



To put the Items together



- **Hardware ToxProtect**
- **Fish/Testorganism**
- **Alarm Software**



Application Measuring Station



Real alarm from pig farm faeces in Latvia



Application Water Quality Monitoring

**Hydropower dam
upside Tehran, Iran**





Summary

ToxProtect is a qualified tool for water quality survey. It detects contaminations in short times.

ToxProtect combines the ability of fish toxicity assessment with a superior alarm recognition.

ToxProtect includes an alarm verification check.

ToxProtect reduces costs and risks.

Thank you

